

SAFETY DATA SHEET

According to Regulation 1907/2006/CE

Doc. N° 401034 rev.1
Date of issue: 19/01/2015

1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1 Product identifier

Product name : **Amies Transport Medium**
Product Number : **401034**

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses : Laboratory chemicals, culture medium for microbiology

1.3 Details of the supplier of the safety data sheet

Company : Biolife Italiana S.r.l.
Viale Monza 272, 20128 Milano Italia
Tel : 0039 02 252091
Fax: 0039 02 2576428
E-mail: mktg@biolifeitaliana.it

1.4 Emergency telephone number

Emergency Phone : 0039 02-6610-1029 (Centro Antiveneni Niguarda Ca'
Granda- Milano)

2. HAZARDS IDENTIFICATION

2.1 Classification of the substance or mixture

Not a hazardous substance or mixture according to
Regulation (EC) No. 1272/2008.
This substance is not classified as dangerous according to
Directive 67/548/EEC.

2.2 Label elements

Labelling according Regulation (EC) No 1272/2008 [CLP]

Pictogram none
Signal word none
Hazard statement(s) none
Precautionary statement(s) none
Supplemental Hazard
Statements none
Caution - this mixture contains a substance not yet fully tested.

2.3 Other hazards – none

3. COMPOSITION/INFORMATION ON INGREDIENTS

3.2 Mixtures

Component	Classification	Concentration
Sodium mercaptoacetate		
CAS-No 367-51-1 CE-No 206-696-4	Acute Tox. 4; Skin Irrit. 2; Eye Irrit. 2; H302, H315, H319 ----- Xn, R22 - R38	< 10 %

For the full text of the H-Statements and R-Phrases mentioned in this Section, see Section 16

4. FIRST AID MEASURES

Biolife Italiana S.r.l., Viale Monza 272, 20128 Milan, Italy. Tel. n°++39 02 25209.1, Fax n° ++39 02 2576428
E-mail: mktg@biolifeitaliana.it; Web Site: www.biolifeitaliana.it

4.1 Description of first aid measures

Consult a physician. Show this safety data sheet to the doctor in attendance.

If inhaled

If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact

Wash off with soap and plenty of water. Consult a physician.

In case of eye contact

Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed

Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

4.2 Most important symptoms and effects, both acute and delayed

To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.

4.3 Indication of any immediate medical attention and special treatment needed

no data available

5. FIREFIGHTING MEASURES**5.1 Extinguishing media****Suitable extinguishing media**

Use water spray, alcohol-resistant foam, dry chemical or carbon dioxide.

5.2 Special hazards arising from the substance or mixture

no data available

5.3 Advice for firefighters

Wear self contained breathing apparatus for fire fighting if necessary.

5.4 Further information

no data available

6. ACCIDENTAL RELEASE MEASURES**6.1 Personal precautions, protective equipment and emergency procedures**

Use personal protective equipment. Avoid dust formation. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Avoid breathing dust.

6.2 Environmental precautions

Do not let product enter drains.

6.3 Methods and materials for containment and cleaning up

Pick up and arrange disposal without creating dust. Sweep up and shovel. Keep in suitable, closed containers for disposal.

6.4 Reference to other sections

For disposal see section 13.

7. HANDLING AND STORAGE**7.1 Precautions for safe handling**

Avoid contact with skin and eyes. Avoid formation of dust and aerosols.

Provide appropriate exhaust ventilation at places where dust is formed.

7.2 Conditions for safe storage, including any incompatibilities

Biolife Italiana S.r.l., Viale Monza 272, 20128 Milan, Italy. Tel. n° ++39 02 25209.1, Fax n° ++39 02 2576428

E-mail: mktg@biolifeitaliana.it; Web Site: www.biolifeitaliana.it

Store in cool place. Keep container tightly closed in a dry and well-ventilated place.

7.3 Specific end uses

Apart from the uses mentioned in section 1.2.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1 Control parameters

Components with workplace control parameters

It doesn't contain substances with occupational exposure limit value.

8.2 Exposure controls

Appropriate engineering controls

Handle in accordance with good industrial hygiene and safety practice. Wash hands before breaks and at the end of workday.

Personal protective equipment

Eye/face protection

Safety glasses with side-shields conforming to EN166 Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

Skin protection

Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices. Wash and dry hands. The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Body Protection

impervious clothing, The type of protective equipment must be selected according to the concentration and amount of the dangerous substance at the specific workplace.

Respiratory protection

For nuisance exposures use type P95 (US) or type P1 (EU EN 143) particle respirator. For higher level protection use type OV/AG/P99 (US) or type ABEK-P2 (EU EN 143) respirator cartridges. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

9. PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

a) Appearance	powder black
b) Odour	no data available
c) Odour Threshold	no data available
d) pH	6.9 – 7.3 at 25°C
e) Melting point/freezing point	no data available
f) Initial boiling point and boiling range	no data available
g) Flash point	no data available
h) Evaporation rate	no data available
i) Flammability (solid, gas)	no data available
j) Upper/lower	

flammability or explosive limits	no data available
k) Vapour pressure	no data available
l) Vapour density	no data available
m) Relative density	no data available
n) Water solubility	no data available
o) Partition coefficient: noctanol/water	no data available
p) Auto-ignition temperature	no data available
q) Decomposition temperature	no data available
r) Viscosity	no data available
s) Explosive properties	no data available
t) Oxidizing properties	no data available

9.2 Other safety information no data available

10. STABILITY AND REACTIVITY

10.1 Reactivity	no data available
10.2 Chemical stability	no data available
10.3 Possibility of hazardous reactions	no data available
10.4 Conditions to avoid	no data available
10.5 Incompatible materials	no data available
10.6 Hazardous decomposition products	no data available

11. TOXICOLOGICAL INFORMATION

11.1 Information on toxicological effects

Acute toxicity	no data available
Skin corrosion/irritation	no data available
Serious eye damage/eye irritation	no data available
Respiratory or skin sensitization	no data available
Germ cell mutagenicity	no data available
Carcinogenicity	IARC: No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
Reproductive toxicity	no data available
Specific target organ toxicity - single exposure	no data available
Specific target organ toxicity - repeated exposure	no data available
Potential health effects	
Inhalation	May be harmful if inhaled. Causes respiratory tract irritation.

Ingestion	Harmful if swallowed.
Skin	May be harmful if absorbed through skin. Causes skin irritation.
Eyes	Causes serious eye irritation.
Signs and Symptoms of Exposure	To the best of our knowledge, the chemical, physical, and toxicological properties have not been thoroughly investigated.
Additional Information	RTECS: Not available

12. ECOLOGICAL INFORMATION

12.1 Toxicity	no data available
12.2 Persistence and degradability	no data available
12.3 Bioaccumulative potential	no data available
12.4 Mobility in soil	no data available
12.5 Results of PBT and vPvB assessment	no data available
12.6 Other adverse effects	no data available

13. DISPOSAL CONSIDERATIONS
13.1 Waste treatment methods

Product Offer surplus and non-recyclable solutions to a licensed disposal company. Contact a licensed professional waste disposal service to dispose of this material. Dissolve or mix the material with a combustible solvent and burn in a chemical incinerator equipped with an afterburner and scrubber.

Contaminated packaging Dispose of as unused product.

14. TRANSPORT INFORMATION

14.1 UN number	ADR/RID: - IMDG: - IATA: -
14.2 UN proper shipping name	ADR/RID: Not dangerous goods IMDG: Not dangerous goods IATA: Not dangerous goods
14.3 Transport hazard class(es)	ADR/RID: - IMDG: - IATA: -
14.4 Packaging group	ADR/RID: - IMDG: - IATA: -
14.5 Environmental hazards	ADR/RID: no IMDG Marine pollutant: no IATA: no
14.6 Special precautions for user	no data available

15. REGULATORY INFORMATION

This safety datasheet complies with the requirements of Regulation (EC) No. 1907/2006.

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Biolife Italiana S.r.l., Viale Monza 272, 20128 Milan, Italy. Tel. n° ++39 02 25209.1, Fax n° ++39 02 2576428
E-mail: mktg@biolifeitaliana.it; Web Site: www.biolifeitaliana.it

no data available

15.2 Chemical Safety Assessment

no data available

16. OTHER INFORMATION**Text of H-code(s) and R-phrases mentioned in Section 3**

Acute Tox.	Acute toxicity
Eye Irrit.	Eye irritation
H302	Harmful if swallowed.
H315	Causes skin irritation.
H319	Causes serious eye irritation.
Skin Irrit.	Skin irritation
Xn	Harmful
R22	Harmful if swallowed.
R38	Irritating to skin.

Further information**Restrictions for use**

No data available

Training advice

No data available

References

No data available

The information in this MSDS was obtained from current and reliable sources. Anyone using this information is solely responsible for the accuracy and applicability of this information to a particular use or situation. Biolife Italiana S.r.l. does not in any way warrant or imply the applicability, viability or use of this information to any person or for use in any situation.